

Datasheet: VMA00011KT

Description:	CYTOKERATIN 18/19 ANTIBODY WITH CONTROL LYSATE
Specificity:	CYTOKERATIN 18/19
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Isotype:	HuCAL Fab bivalent
Quantity:	2 Westerns

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	A bivalent human recombinant Fab selected from the HuCAL [®] GOLD phage display library, expressed in <i>E. coli</i> . This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This product is supplied as a liquid.
Preparation	20µl Recombinant human monoclonal antibody purified by metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Recombinant human cytokeratin 18

**External Database
Links**

UniProt:

[P05783](#)

[Related reagents](#)

[P08727](#)

[Related reagents](#)

Entrez Gene:

[3875](#)

KRT18

[Related reagents](#)

[3880](#)

KRT19

[Related reagents](#)

Synonyms

CYK18

Specificity

Human anti Human cytokeratin 18/19 antibody recognizes human cytokeratin 18, a amino acid type I intermediate filament protein expressed by epithelial cells. The intermediate filaments are expressed as heterodimers in epithelial tissues. Cytokeratin 18 (CK18) is often co-expressed with cytokeratin 8 (CK8). The CK18/CK8 heterodimer is the major cytokeratin pair found in simple epithelia such as that found in the liver, intestine and breast ducts. CK18 plays an important role in maintaining epithelial cell integrity.

This clone has also been demonstrated to cross-react with cytokeratin 19. Human anti Human cytokeratin 18/19 antibody recognizes cytokeratin 18 as a band of ~48 kDa and cytokeratin 19 as a band of ~44 kDa in western blotting under reducing conditions.

Western Blotting

Anti cytokeratin 18/19 antibody detects bands of approximately 44 kDa in HeLa cell lysates.

Instructions For Use

Please refer to the [PrecisionAb western blotting protocol](#). For additional information on secondary antibody dilution and exposure time see product web page.

Lysate Composition

400µg HeLa lysate lyophilized in RIPA buffer.

Lysate Reconstitution

- If using DDT reconstitute the lyophilized lysate with 190µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.
- If using BME reconstitute the lyophilized lysate with 180µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.
Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.

Storage

Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.

Guarantee

As supplied, 12 months from date of despatch.

Acknowledgements

PrecisionAb™ is a trademark of Bio-Rad Laboratories. His-tag is a registered trademark of EMD Biosciences.

Health And Safety

Material Safety Datasheet documentation #10040 #10561 available at:

Information Antibody (10040): <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>
Lysate Material (10561): <https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Goat Anti Human IgG F(ab')₂ (STAR126...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M351019:190315'

Printed on 09 Feb 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)