

Datasheet: VMA00083

Description:	MOUSE ANTI HIF2 ALPHA
Specificity:	HIF2 ALPHA
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Mouse monoclonal antibody purified by affinity chromatography on Protein G from tissue culture supernatant.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)
1% Bovine Serum Albumin
<50% Glycerol.

Immunogen	Recombinant protein corresponding to an internal region of human HIF2 alpha
External Database Links	<p>UniProt: Q99814 Related reagents</p> <p>Entrez Gene: 2034 EPAS1 Related reagents</p>
Synonyms	BHLHE73, HIF2A, MOP2, PASD2
Specificity	<p>Mouse anti Human HIF2 alpha antibody recognizes HIF2 alpha, also known as EPAS-1, basic-helix-loop-helix-PAS protein MOP2, class E basic helix-loop-helix protein 73, bHLHe73, HIF-1-alpha-like factor, HLF, hypoxia-inducible factor 2-alpha, HIF-2-alpha, HIF2-alpha, member of PAS protein 2 and PAS domain-containing protein 2.</p> <p>HIF 2 alpha is a transcription factor involved in the induction of genes regulated by oxygen, which is induced as oxygen levels fall. Encoded by the EPAS1 gene, HIF2 alpha contains a basic-helix-loop-helix domain protein dimerization domain as well as a domain found in proteins in signal transduction pathways which respond to oxygen levels. Mutations in HIF2 alpha are associated with erythrocytosis familial type 4 (provided by RefSeq, Nov 2009).</p> <p>Mouse anti Human HIF2 alpha antibody detects a band of 117 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti HIF2 alpha detects a band of approximately 117 kDa in HeLa cell lysates.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10048 available at: https://www.bio-rad-antibodies.com/SDS/VMA00083 Antibody (10048)
Regulatory	For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

To
find a

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M397213:220616'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)