

## Datasheet: VMA00083

<b>Description:</b>	MOUSE ANTI HIF2 ALPHA
<b>Specificity:</b>	HIF2 ALPHA
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse, Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Preparation

Mouse monoclonal antibody purified by affinity chromatography on Protein G from tissue culture supernatant

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)  
1% Bovine Serum Albumin  
<50% Glycerol

<b>Immunogen</b>	Recombinant protein corresponding to an internal region of human HIF2 alpha
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">Q99814</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">2034</a>    EPAS1    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	BHLHE73, HIF2A, MOP2, PASD2
<b>Specificity</b>	<p><b>Mouse anti Human HIF2 alpha antibody</b> recognizes HIF2 alpha, also known as EPAS-1, basic-helix-loop-helix-PAS protein MOP2, class E basic helix-loop-helix protein 73, bHLHe73, HIF-1-alpha-like factor, HLF, hypoxia-inducible factor 2-alpha, HIF-2-alpha, HIF2-alpha, member of PAS protein 2 and PAS domain-containing protein 2.</p> <p>HIF 2 alpha is a transcription factor involved in the induction of genes regulated by oxygen, which is induced as oxygen levels fall. Encoded by the EPAS1 gene, HIF2 alpha contains a basic-helix-loop-helix domain protein dimerization domain as well as a domain found in proteins in signal transduction pathways which respond to oxygen levels. Mutations in HIF2 alpha are associated with erythrocytosis familial type 4 (provided by RefSeq, Nov 2009).</p> <p>Mouse anti Human HIF2 alpha antibody detects a band of 117 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti HIF2 alpha detects a band of approximately 117 kDa in HeLa cell lysates.
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10048 available at: Antibody (10048): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

'M374218:201023'

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

**Printed on 29 Oct 2020**

---

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)