

Datasheet: VMA00117

BATCH NUMBER 080615

Description:	MOUSE ANTI CYTOCHROME P450 3A4
Specificity:	CYTOCHROME P450 3A4
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody prepared by affinity chromatography on Protein G from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Purified recombinant fragment of human cytochrome P450 3A4 expressed in <i>E. coli</i>
External Database Links	UniProt: P08684 Related reagents

Entrez Gene:

[1576](#) CYP3A4 [Related reagents](#)

Synonyms	CYP3A3
-----------------	--------

Specificity	<p>Mouse anti Human cytochrome P450 3A4 antibody recognizes cytochrome P450 3A4, also known as 1,8-cineole 2-exo-monooxygenase, P450-III, steroid inducible, albendazole monooxygenase, albendazole sulfoxidase, cytochrome P450 3A3, cytochrome P450 HLP, cytochrome P450 NF-25, cytochrome P450, subfamily IIIA (niphedipine oxidase), polypeptide 3, cytochrome P450, subfamily IIIA (niphedipine oxidase), polypeptide 4, cytochrome P450-PCN1, glucocorticoid-inducible P450, nifedipine oxidase, quinine 3-monooxygenase and taurochenodeoxycholate 6-alpha-hydroxylase.</p> <p>Encoded by the CYP3A4 gene, cytochrome P450 3A4 is a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases that catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and its expression is induced by glucocorticoids and some pharmacological agents. This enzyme is involved in the metabolism of approximately half the drugs in use today, including acetaminophen, codeine, cyclosporin A, diazepam and erythromycin. The enzyme also metabolizes some steroids and carcinogens. CYP3A4 is part of a cluster of cytochrome P450 genes on chromosome 7q21.1. Previously another CYP3A gene, CYP3A3, was thought to exist; however, it is now thought that this sequence represents a transcript variant of CYP3A4. Alternatively spliced transcript variants encoding different isoforms have been identified (provided by RefSeq, Feb 2011).</p> <p>Mouse anti Human cytochrome P450 3A4 antibody detects a band of 50 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
--------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Western Blotting	Anti cytochrome P450 3A4 detects a band of approximately 50 kDa in human liver tissue cell lysates
-------------------------	----------------------------------------------------------------------------------------------------

Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
----------------	-----------------------------------------------------------------

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
-------------------------	-----------------------------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/VMA00117 Antibody (10040)
--------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M369915:200529'

Printed on 10 Apr 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)