

Datasheet: VMA00132

BATCH NUMBER 240615

Description:	MOUSE ANTI p130Cas
Specificity:	p130Cas
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	OTI3A11
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Mouse monoclonal antibody purified by affinity chromatography from ascites.

Buffer Solution

Phosphate buffered saline

Preservative

0.09% Sodium Azide (NaN₃)

Stabilisers	1% Bovine Serum Albumin 50% Glycerol
Immunogen	Full length recombinant protein of human p130Cas (NP_055382) produced in HEK293T cells
External Database Links	<p>UniProt: P56945 Related reagents</p> <p>Entrez Gene: 9564 BCAR1 Related reagents</p>
Synonyms	CAS, CASS1, CRKAS
Specificity	<p>Mouse anti Human p130Cas antibody, clone OTI3A11 recognizes p130Cas, also known as Cas scaffolding protein family member 1, Crk-associated substrate p130Cas and breast cancer anti-estrogen resistance protein 1.</p> <p>p130Cas, encode by the BCAR1 gene, is an Src family kinase substrate involved in various cellular events, including migration, survival, transformation, and invasion (Sawada <i>et al.</i> 2006).</p> <p>Mouse anti Human p130Cas antibody detects a band of 130 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti p130Cas detects a band of approximately 130 kDa in MCF7 cell lysates.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10048 available at: https://www.bio-rad-antibodies.com/SDS/VMA00132 Antibody (10048)
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M369928:200529'

Printed on 13 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)