

Datasheet: VMA00142

Description:	MOUSE ANTI PSMB4
Specificity:	PSMB4
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	OTI5B5
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody purified by affinity chromatography from tissue culture supernatant
Buffer Solution	Phosphate buffered saline.
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 25% Glycerol
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Full length recombinant protein of human PSMB4 (NP_002787) produced in HEK293T cell

**External Database
Links**

UniProt:

[P28070](#)

[Related reagents](#)

Entrez Gene:

[5692](#)

PSMB4

[Related reagents](#)

Synonyms

PROS26

Specificity

Mouse anti Human PSMB4 antibody, clone OTI5B5 recognizes PSMB4, also known as 26 kDa prosomal protein, hsBPROS26, macropain beta chain, multicatalytic endopeptidase complex beta chain, proteasome beta chain, proteasome chain 3, proteasome subunit HsN3 and proteasome subunit beta type-4.

The proteasome is a multicatalytic proteinase complex with a highly ordered ring-shaped 20S core structure. The core structure is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. PSMB4 encodes a member of the proteasome B-type family, also known as the T1B family, that is a 20S core beta subunit (provided by RefSeq, Jul 2008).

Mouse anti Human PSMB4 antibody detects a band of 28 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti PSMB4 detects a band of approximately 28 kDa in Jurkat cell lysates.

Storage

This product is shipped at ambient temperature.
Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee

12 months from date of despatch.

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories.

**Health And Safety
Information**

Material Safety Datasheet documentation #10048 available at:
<https://www.bio-rad-antibodies.com/SDS/VMA00142>
Antibody (10048)

Regulatory

For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M441712:250523'

Printed on 23 May 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)