

Datasheet: VMA00230

BATCH NUMBER 160321

Description:	MOUSE ANTI HEXOKINASE 1
Specificity:	HEXOKINASE 1
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	7A7
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody prepared by affinity chromatography on Protein G from ascites
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Purified recombinant fragment of human hexokinase 1 expressed in <i>E. coli</i>
External Database Links	<p>UniProt: P19367 Related reagents</p> <p>Entrez Gene: 3098 HK1 Related reagents</p>
Specificity	<p>Mouse anti Human hexokinase 1 antibody recognizes hexokinase 1, also known as HK I, brain form hexokinase, glycolytic enzyme and hexokinase type I.</p> <p>Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. Encoded by the HK1 gene, hexokinase 1 represents a ubiquitous form of hexokinase which localizes to the outer membrane of mitochondria. Mutations in HK1 have been associated with hemolytic anemia due to hexokinase deficiency. Alternative splicing of HK1 results in five transcript variants which encode different isoforms, some of which are tissue-specific. Each isoform has a distinct N-terminus; the remainder of the protein is identical among all the isoforms. A sixth transcript variant has been described, but due to the presence of several stop codons, it is not thought to encode a protein (provided by RefSeq, Apr 2009).</p> <p>Mouse anti Human hexokinase 1 antibody detects a band of 102 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti hexokinase 1 detects a band of approximately 102 kDa in K562 cell lysates
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/VMA00230
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M380374:210511'

Printed on 09 Jul 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)