

Datasheet: VMA00242

**BATCH NUMBER 170976**

<b>Description:</b>	MOUSE ANTI STAT3
<b>Specificity:</b>	STAT3
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	3B5
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Preparation

Mouse monoclonal antibody prepared by affinity chromatography on Protein G from ascites.

### Buffer Solution

Phosphate buffered saline.

### Preservative

0.09% Sodium Azide (NaN<sub>3</sub>).

## Stabilisers

---

**Approx. Protein Concentrations** IgG concentration 0.5 mg/ml.

---

**Immunogen** Purified recombinant fragment of human STAT3 expressed in *E. coli*

---

## External Database Links

**UniProt:**

[P40763](#)    [Related reagents](#)

**Entrez Gene:**

[6774](#) STAT3    [Related reagents](#)

---

**Synonyms** APRF

---

## Specificity

**Mouse anti Human STAT3 antibody** recognizes signal transducer and activator of transcription 3, also known as DNA-binding protein APRF, acute-phase response factor and signal transducer and activator of transcription 3.

Encoded by STAT3 gene, signal transducer and activator of transcription 3 is a member of the STAT protein family. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators. This protein is activated through phosphorylation in response to various cytokines and growth factors including IFNs, EGF, IL5, IL6, HGF, LIF and BMP2. This protein mediates the expression of a variety of genes in response to cell stimuli, and thus plays a key role in many cellular processes such as cell growth and apoptosis. The small GTPase Rac1 has been shown to bind and regulate the activity of this protein. PIAS3 protein is a specific inhibitor of this protein. Three alternatively spliced transcript variants encoding distinct isoforms have been described (provided by RefSeq, Jul 2008).

Mouse anti Human STAT3 antibody detects a band of 84 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting** Anti STAT3 detects a band of approximately 84 kDa in LNCaP cell lysates.

---

**Storage** Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee** 12 months from date of despatch.

---

**Acknowledgements** PrecisionAb is a trademark of Bio-Rad Laboratories.

---

**Health And Safety Information** Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/VMA00242>

---

**Regulatory** For research purposes only.

---

## Related Products

## Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Product inquiries: [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M429458:240410'

Printed on 09 Jul 2025

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)