

Datasheet: VMA00245

Description:	MOUSE ANTI DNA METHYLTRANSFERASE 1
Specificity:	DNA METHYLTRANSFERASE 1
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody prepared by affinity chromatography on Protein G from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Recombinant protein corresponding to amino acids 1448-1594 of human DNA methyltransferase 1 expressed in <i>E. coli</i>
External Database Links	UniProt: P26358 Related reagents

Entrez Gene:[1786](#) DNMT1 [Related reagents](#)

Synonyms	AIM, CXXC9, DNMT
-----------------	------------------

Specificity	<p>Mouse anti Human DNA methyltransferase 1 antibody recognizes the DNA (cytosine-5)-methyltransferase 1 , also known as Dnmt1, CXXC-type zinc finger protein 9, DNA methyltransferase Hsa1, DNA MTase Hsa1 or M.Hsa1 MCMT.</p> <p>DNA (cytosine-5-)-methyltransferase 1 has a role in the establishment and regulation of tissue-specific patterns of methylated cytosine residues. Aberrant methylation patterns are associated with certain human tumors and developmental abnormalities. Two transcript variants encoding different isoforms have been found for DNMT1 (provided by RefSeq, Aug 2008).</p> <p>Mouse anti Human DNA methyltransferase 1 antibody detects a band of 180 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
--------------------	--

Western Blotting	Anti DNA methyltransferase 1 detects a band of approximately 205 kDa in Jurkat cell lysates
-------------------------	---

Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
----------------	---

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
-------------------------	---

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M377135:210212'

Printed on 02 Mar 2021

