

Datasheet: VMA00247

**BATCH NUMBER 151112**

<b>Description:</b>	MOUSE ANTI E3 UBIQUITIN-PROTEIN LIGASE CBL
<b>Specificity:</b>	E3 UBIQUITIN-PROTEIN LIGASE CBL
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	3B12
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunoprecipitation	▪			
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Preparation

Mouse monoclonal antibody prepared by affinity chromatography on Protein G from ascites

### Buffer Solution

Phosphate buffered saline

<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	Purified recombinant fragment of human E3 ubiquitin-protein ligase CBL expressed in <i>E. coli</i>
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P22681</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">867</a> CBL    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	CBL2, RNF55
<b>Specificity</b>	<p><b>Mouse anti Human E3 ubiquitin-protein ligase CBL antibody</b> recognizes the E3 ubiquitin-protein ligase CBL, also known as Cas-Br-M (murine) ecotropic retroviral transforming sequence, E3 ubiquitin-protein ligase CBL, RING finger protein 55, casitas B-lineage lymphoma proto-oncogene, fragile site, folic acid type, rare, fra(11)(q23.3), oncogene CBL2, proto-oncogene c-Cbl and signal transduction protein CBL.</p> <p>Encoded by the CBL gene, E3 ubiquitin-protein ligase CBL is a proto-oncogene that encodes a RING finger E3 ubiquitin ligase. E3 ubiquitin-protein ligase CBL is one of the enzymes required for targeting substrates for degradation by the proteasome. This protein mediates the transfer of ubiquitin from ubiquitin conjugating enzymes (E2) to specific substrates. This protein also contains an N-terminal phosphotyrosine binding domain that allows it to interact with numerous tyrosine-phosphorylated substrates and target them for proteasome degradation. As such it functions as a negative regulator of many signal transduction pathways. CBL has been found to be mutated or translocated in many cancers including acute myeloid leukaemia. Mutations in CBL are also the cause of Noonan syndrome-like disorder (provided by RefSeq, Mar 2012).</p> <p>Mouse anti Human E3 ubiquitin-protein ligase CBL antibody detects a band of 120 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti E3 ubiquitin-protein ligase CBL detects a band of approximately 120 kDa in HEK293 cell lysates
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/VMA00247">https://www.bio-rad-antibodies.com/SDS/VMA00247</a>

## Antibody (10040)

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M393963:220120'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)