

## Datasheet: VMA00254

<b>Description:</b>	MOUSE ANTI ALDEHYDE DEHYDROGENASE 3
<b>Specificity:</b>	ALDEHYDE DEHYDROGENASE 3
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	OT11B6
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Mouse monoclonal antibody purified by affinity chromatography from ascites
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Stabilisers</b>	1% Bovine Serum Albumin 50% Glycerol
<b>Immunogen</b>	Full length recombinant human aldehyde dehydrogenase 3 (NP_000682) produced in HEK293T cells
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P30838</a>   <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">218</a>   ALDH3A1   <a href="#">Related reagents</a></p>

**Synonyms** ALDH3

**Specificity** **Mouse anti Human aldehyde dehydrogenase 3 antibody** recognizes the aldehyde dehydrogenase 3, also known as aldehyde dehydrogenase, dimeric NADP-preferring, aldehyde dehydrogenase isozyme 3, aldehyde dehydrogenase type III, aldehyde dehydrogenase, dimeric NADP-preferring and stomach aldehyde dehydrogenase.

Aldehyde dehydrogenases oxidize various aldehydes to the corresponding acids. They are involved in the detoxification of alcohol-derived acetaldehyde and in the metabolism of corticosteroids, biogenic amines, neurotransmitters, and lipid peroxidation. The enzyme encoded by ALDH3A1 gene forms a cytoplasmic homodimer that preferentially oxidizes aromatic and medium-chain (6 carbons or more) saturated and unsaturated aldehyde substrates. It is thought to promote resistance to UV and 4-hydroxy-2-nonenal-induced oxidative damage in the cornea. The gene is located within the Smith-Magenis syndrome region on chromosome 17. Multiple alternatively spliced variants, encoding the same protein, have been identified (provided by RefSeq, Sep 2008).

Mouse anti Human aldehyde dehydrogenase 3 antibody detects a band of 50 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

**Western Blotting** Anti aldehyde dehydrogenase 3 detects a band of approximately 50 kDa in A549 cell lysates

**Storage** Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

**Guarantee** 12 months from date of despatch

**Acknowledgements** PrecisionAb is a trademark of Bio-Rad Laboratories.

**Health And Safety Information** Material Safety Datasheet documentation #10048 available at: Antibody (10048): <https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf>

**Regulatory** For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

### Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M370031:200529'

Printed on 11 Aug 2020