

## Datasheet: VMA00274

<b>Description:</b>	MOUSE ANTI DNAJA2
<b>Specificity:</b>	DNAJA2
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	OTI3A11
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse, Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Preparation

Mouse monoclonal antibody purified by affinity chromatography from ascites.

### Buffer Solution

Phosphate buffered saline.

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)  
1% Bovine Serum Albumin  
50% Glycerol

<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	Full length recombinant human DNAJA2 (NP_005871) produced in HEK293T cells
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">O60884</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">10294</a> DNAJA2    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	CPR3, HIRIP4
<b>Specificity</b>	<p><b>Mouse anti Human DNAJA2 antibody, clone OT13A11</b> recognizes the DnaJ homolog subfamily A member 2, also known as HIRA interacting protein 4, HIRA-interacting protein 4, cell cycle progression 3 protein, cell cycle progression restoration gene 3 protein, dnaJ homolog subfamily A member 2 and renal carcinoma antigen NY-REN-14.</p> <p>Encoded by DNAJA2 gene, DnaJ homolog subfamily A member 2 belongs to the evolutionarily conserved DNAJ/HSP40 family of proteins, which regulate molecular chaperone activity by stimulating ATPase activity. DNAJ proteins may have up to 3 distinct domains: a conserved 70-amino acid J domain, usually at the N terminus; a glycine/phenylalanine (G/F)-rich region; and a cysteine-rich domain containing 4 motifs resembling a zinc finger domain. The product of DNAJA2 works as a cochaperone of Hsp70s in protein folding and mitochondrial protein import in vitro (provided by RefSeq, Jul 2008).</p> <p>Mouse anti Human DNAJA2 antibody, clone OT13A11 detects a band of 46 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti DNAJA2 detects a band of approximately 46 kDa in HEK293 cell lysates.
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
<b>Guarantee</b>	12 months from date of despatch.
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10048 available at: <a href="https://www.bio-rad-antibodies.com/SDS/VMA00274">https://www.bio-rad-antibodies.com/SDS/VMA00274</a> Antibody (10048)
<b>Regulatory</b>	For research purposes only.

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

## Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M415577:230111'

Printed on 13 Aug 2023

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)