

Datasheet: VMA00282

BATCH NUMBER 151112

Description:	MOUSE ANTI DREBRIN-LIKE PROTEIN
Specificity:	DREBRIN-LIKE PROTEIN
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	OTI7B7
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody purified by affinity chromatography from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 50% Glycerol
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml

Immunogen Full length recombinant human drebrin-like protein (NP_054782) produced in HEK293T cells

External Database Links

UniProt:

[Q9UJU6](#) [Related reagents](#)

Entrez Gene:

[28988](#) DBNL [Related reagents](#)

Synonyms

CMAP, SH3P7

Specificity

Mouse anti Human drebrin-like protein antibody, clone OTI7B7 recognizes drebrin-like protein, also known as HPK1-interacting protein of 55 kDa, SH3 domain-containing protein 7, cervical SH3P7, cervical mucin-associated protein, drebrin-F, drebrin-like protein and src homology 3 domain-containing protein HIP-55.

Mouse anti Human drebrin-like protein antibody, clone OTI7B7 detects a band of 55 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti drebrin-like protein detects a band of approximately 55 kDa in HeLa cell lysates

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee

12 months from date of despatch

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories.

Health And Safety Information

Material Safety Datasheet documentation #10048 available at: <https://www.bio-rad-antibodies.com/SDS/VMA00282>
Antibody (10048)

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M393251:211215'

Printed on 25 Mar 2023
