

Datasheet: VMA00314

BATCH NUMBER 151112

Description:	MOUSE ANTI CYTOKERATIN 18
Specificity:	CYTOKERATIN 18
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Isotype:	IgG2b
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody prepared by affinity chromatography on Protein G from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Recombinant protein corresponding to amino acid 391-483 of human cytokeratin 18 expressed in <i>E. coli</i>
External Database Links	UniProt:

Entrez Gene:

[3875](#) KRT18 [Related reagents](#)

Synonyms CYK18

Specificity **Mouse anti Human cytokeratin 18 antibody** recognizes cytokeratin-18, also known as keratin, type I cytoskeletal 18, CK-18, cell proliferation-inducing gene 46 protein, cell proliferation-inducing protein 46, keratin, type I cytoskeletal 18 and keratin-18.

Encoded by the KRT18 gene cytokeratin 18, together with its filament partner cytokeratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in KRT18 have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for KRT18 (provided by RefSeq, Jul 2008).

Mouse anti Human cytokeratin 18 antibody detects a band of 46 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting Anti cytokeratin 18 detects a band of approximately 46 kDa in MCF7 cell lysates

Storage Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee 12 months from date of despatch

Acknowledgements PrecisionAb is a trademark of Bio-Rad Laboratories.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/VMA00314>
Antibody (10040)

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M370086:200529'

Printed on 25 Mar 2023