

## Datasheet: VMA00319

<b>Description:</b>	MOUSE ANTI LIN28A
<b>Specificity:</b>	LIN28A
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb™ Monoclonal
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**PrecisionAb antibodies have been extensively [validated for the western blot application](#).** The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Mouse monoclonal antibody prepared by affinity chromatography on Protein G from ascites
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Immunogen</b>	Recombinant protein corresponding to amino acid 93-209 of LIN28A expressed in <i>E. coli</i>
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">Q9H9Z2</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">79727</a>    LIN28A    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	CSDD1, LIN28, ZCCHC1
<b>Specificity</b>	<b>Mouse anti Human LIN28A antibody</b> recognizes the protein lin-28 homolog A, also known as RNA-binding protein LIN-28, protein lin-28 homolog A, zinc finger CCHC domain-containing protein 1, zinc finger and CCHC domain containing 1.

Mouse anti Human LIN28A antibody detects a band of 28 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

<b>Western Blotting</b>	Anti LIN28A detects a band of approximately 28 kDa in NTERA-2 cell lysates
-------------------------	--

---

<b>Instructions For Use</b>	Please refer to the <a href="#">PrecisionAb western blotting protocol</a> . For additional information on secondary antibody dilution and exposure time see product web page.
-----------------------------	---

---

<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
----------------	---

---

<b>Guarantee</b>	As supplied, 12 months from date of despatch.
------------------	---

---

<b>Acknowledgements</b>	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
-------------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: Antibody (10040): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</a>
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M337563:181217'

**Printed on 15 Mar 2019**