

Datasheet: VMA00345

BATCH NUMBER 151215

Description:	MOUSE ANTI NUCLEOSIDE DIPHOSPHATE KINASE A
Specificity:	NUCLEOSIDE DIPHOSPHATE KINASE A
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Isotype:	IgG2b
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Mouse monoclonal antibody prepared by affinity chromatography on Protein G from tissue culture supernatant.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃).

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant human nucleoside diphosphate kinase A
External Database Links	<p>UniProt: P15531 Related reagents</p> <p>Entrez Gene: 4830 NME1 Related reagents</p>
Synonyms	NDPKA, NM23
Specificity	<p>Mouse anti Human nucleoside diphosphate kinase A antibody recognizes nucleoside diphosphate kinase A, also known as NDP kinase A, granzyme A-activated DNase, metastasis inhibition factor nm23, non-metastatic cells 1 and tumor metastatic process-associated protein.</p> <p>The NME1 gene was identified because of its reduced mRNA transcript levels in highly metastatic cells. Nucleoside diphosphate kinase (NDK) exists as a hexamer composed of 'A' (encoded by NME1) and 'B' (encoded by NME2) isoforms. Mutations in NME1 have been identified in aggressive neuroblastomas. Two transcript variants encoding different isoforms have been found for NME1. Co-transcription of NME1 and the neighboring downstream gene (NME2) generates naturally-occurring transcripts (NME1-NME2), which encodes a fusion protein comprised of sequence sharing identity with each individual gene product (provided by RefSeq, Jul 2008).</p> <p>Mouse anti Human nucleoside diphosphate kinase A antibody detects a band of 17 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti nucleoside diphosphate kinase A detects a band of approximately 17 kDa in HeLa cell lysates.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch.
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/VMA00345 Antibody (10040)
Regulatory	For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL \(MCA691\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M404440:220825'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)