

Datasheet: VMA00379

Description:	MOUSE ANTI HNF4A
Specificity:	HNF4A
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	OTI2H2
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody purified by affinity chromatography from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 50% Glycerol
Immunogen	Recombinant protein fragment corresponding to aa 1-315 of human HNF4A (NP_849180) produced in <i>E.coli</i>
External Database	UniProt:

Links

[P41235](#) [Related reagents](#)

Entrez Gene:

[3172](#) HNF4A [Related reagents](#)

Synonyms

HNF4, NR2A1, TCF14

Specificity

Mouse anti Human HNF4A antibody recognizes hepatocyte nuclear factor 4-alpha, also known as TCF-14, nuclear receptor subfamily 2 group A member 1 and transcription factor 14.

Encoded by the HNF4A gene, hepatocyte nuclear factor 4-alpha is a nuclear transcription factor which binds DNA as a homodimer. Hepatocyte nuclear factor 4-alpha controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. HNF4A may play a role in development of the liver, kidney, and intestines. Mutations in HNF4A have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of HNF4A results in multiple transcript variants encoding several different isoforms (provided by RefSeq, Apr 2012).

Mouse anti Human HNF4A antibody detects a band of 52 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti HNF4A detects a band of approximately 52 kDa in HepG2 cell lysates

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee

12 months from date of despatch

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories.

Health And Safety Information

Material Safety Datasheet documentation #10048 available at: Antibody (10048): <https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M370145:200529'

