

Datasheet: VMA00384KT

Description:	VISA ANTIBODY WITH CONTROL LYSATE
Specificity:	VISA
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	OTI8A11
Isotype:	IgG1
Quantity:	2 Westerns

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	20µl Mouse monoclonal antibody purified by affinity chromatography from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 50% Glycerol
Immunogen	Recombinant protein fragment corresponding to aa 1-247 of human VISA (NP_065797) produced in <i>E.coli</i>
External Database	UniProt:

Links

[Q7Z434](#) [Related reagents](#)

Entrez Gene:

[57506](#) MAVS [Related reagents](#)

Synonyms

IPS1, KIAA1271, VISA

Specificity

Mouse anti Human VISA antibody recognizes virus-induced signaling adaptor (VISA), also known as MAVS, CARD adapter inducing interferon beta, interferon beta promoter stimulator protein 1, mitochondrial antiviral-signaling protein and putative NF-kappa-B-activating protein 031N.

Encoded by the MAVS gene, VISA is an intermediary protein necessary in the virus-triggered beta interferon signaling pathways. It is required for activation of transcription factors which regulate expression of beta interferon and contributes to antiviral immunity. Multiple transcript variants encoding different isoforms have been found for MAVS or VISA (provided by RefSeq, Sep 2011).

Mouse anti Human VISA antibody detects a band of 75 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti VISA detects a band of approximately 75 kDa in A431 cell lysates

Instructions For Use

Please refer to the [PrecisionAb western blotting protocol](#). For additional information on secondary antibody dilution and exposure time see product web page.

Lysate Composition

400µg A431 lysate lyophilized in RIPA buffer

Lysate Reconstitution

- If using DDT reconstitute the lyophilized lysate with 190µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.
- If using BME reconstitute the lyophilized lysate with 180µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.
Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.

Storage

Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.

Guarantee

As supplied, 12 months from date of despatch

Acknowledgements

PrecisionAb™ is a trademark of Bio-Rad Laboratories

Health And Safety Information

Material Safety Datasheet documentation #10048 #10561 available at:
Antibody (10048): <https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf>
Lysate Material (10561): <https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf>

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at bio-rad-antibodies.com

'M351629:190318'

Printed on 10 Feb 2021
