

Datasheet: VMA00404

Description:	MOUSE ANTI CREATINE PHOSPHOKINASE
Specificity:	CREATINE PHOSPHOKINASE
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	E05/10G2
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Mouse monoclonal antibody purified by affinity chromatography on Protein G from tissue culture supernatant.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃).

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml.
Immunogen	<i>E.coli</i> -derived recombinant human creatine phosphokinase
External Database Links	<p>UniProt:</p> <p>P12277 Related reagents</p> <p>P06732 Related reagents</p> <p>Entrez Gene:</p> <p>1152 CKB Related reagents</p> <p>1158 CKM Related reagents</p>
Synonyms	CKBB, CKMM
Specificity	<p>Mouse anti Human creatine phosphokinase antibody, clone E05/10G2 recognizes creatine phosphokinase which reversibly catalyzes the transfer of phosphate between ATP and creatine to produce creatine phosphate. In cells, the cytosolic creatine phosphokinase enzymes consist of two subunits, which can be either B (brain type) or M (muscle type). There are three different isoenzymes with homodimers or heterodimers of the M and B types: CKMM or CKBB and CKMB respectively. CKBB is found in the brain while CKMB occurs in heart tissue (provided by RefSeq, Jul 2008).</p> <p>Mouse anti Human creatine phosphokinase antibody, clone E05/10G2 detects a band of 42 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti creatine phosphokinase detects a band of approximately 42 kDa in HEK293 cell lysates.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch.
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/VMA00404 Antibody (10040)
Regulatory	For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M417497:230331'

Printed on 13 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)