

Datasheet: VMA00426

Description:	MOUSE ANTI MGMT
Specificity:	MGMT
Format:	Purified
Product Type:	PrecisionAb™ Monoclonal
Clone:	UMAB56
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody purified by affinity chromatography from ascites
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide (NaN ₃)
Stabilisers	1% Bovine Serum Albumin 50% Glycerol
Immunogen	Full length recombinant human MGMT (NP_002403) produced in <i>E.coli</i>
External Database Links	<p>UniProt: P16455 Related reagents</p> <p>Entrez Gene: 4255 MGMT Related reagents</p>
Specificity	Mouse anti Human MGMT antibody recognizes methylguanine-DNA methyltransferase (MGMT), also known as 6-O-methylguanine-DNA methyltransferase, O-6-methylguanine-

DNA-alkyltransferase and methylated-DNA--protein-cysteine methyltransferase.

Alkylating agents are potent carcinogens that can result in cell death, mutation and cancer. The protein encoded by the MGMT gene is a DNA repair protein that is involved in cellular defense against mutagenesis and toxicity from alkylating agents. The protein catalyzes transfer of methyl groups from O(6)-alkylguanine and other methylated moieties of the DNA to its own molecule, which repairs the toxic lesions. Methylation of the promoter has been associated with several cancer types, including colorectal cancer, lung cancer, lymphoma and glioblastoma (provided by RefSeq, Sep 2015).

Mouse anti Human MGMT antibody detects a band of 23 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting	Anti MGMT detects a band of approximately 23 kDa in Jurkat cell lysates
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10048 available at: Antibody (10048): https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M337660:181217'

Printed on 02 Jan 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)