

Datasheet: VMA00437KT

Description:	EPLIN ANTIBODY WITH CONTROL LYSATE
Specificity:	EPLIN
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	OTI1A3
Isotype:	IgG1
Quantity:	2 Westerns

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	20µl Mouse monoclonal antibody purified by affinity chromatography from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 50% Glycerol
Immunogen	Recombinant protein fragment corresponding to aa 1-275 of human EPLIN or LIMA1 (NP_057441) produced in <i>E.coli</i>
External Database	UniProt:

Links[Q9UHB6](#) [Related reagents](#)**Entrez Gene:**[51474](#) LIMA1 [Related reagents](#)**Synonyms**

EPLIN, SREBP3

Specificity

Mouse anti Human EPLIN antibody recognizes epithelial protein lost in neoplasm (EPLIN), also known as LIM domain and actin-binding protein 1, LIMA1 and sterol regulatory element binding protein 3.

Encoded by the EPLIN / LIMA1 gene, epithelial protein lost in neoplasm is a cytoskeleton-associated protein that inhibits actin filament depolymerization and cross-links filaments in bundles. It is downregulated in some cancer cell lines. Alternatively spliced transcript variants encoding different isoforms have been described for LIMA1, and expression of some of the variants maybe independently regulated (provided by RefSeq, Aug 2011).

Mouse anti Human EPLIN antibody detects a band of 110 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti EPLIN detects a band of approximately 110 kDa in HEK293 cell lysates

Instructions For Use

Please refer to the [PrecisionAb western blotting protocol](#). For additional information on secondary antibody dilution and exposure time see product web page.

Lysate Composition

400µg HEK293 lysate lyophilized in RIPA buffer

Lysate Reconstitution

- If using DDT reconstitute the lyophilized lysate with 190µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.

- If using BME reconstitute the lyophilized lysate with 180µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.

Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.

Storage

Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.

Guarantee

As supplied, 12 months from date of despatch

Acknowledgements

PrecisionAb™ is a trademark of Bio-Rad Laboratories

Health And Safety Information

Material Safety Datasheet documentation #10048 #10561 available at:

Antibody (10048): <https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf>Lysate Material (10561): <https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf>**Regulatory**

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at bio-rad-antibodies.com

'M351085:190315'

Printed on 10 Feb 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)