

Datasheet: VMA00444

BATCH NUMBER 160226

Description:	HUMAN ANTI CYTOKERATIN 18
Specificity:	CYTOKERATIN 18
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Isotype:	HuCAL Fab bivalent
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Product Form

A bivalent human recombinant Fab selected from the HuCAL® phage display library, expressed in *E. coli*. This Fab fragment is bivalent by dimerization of the bacterial alkaline phosphatase (BAP) fusion protein. The BAP enzymatic activity is silenced by mutation. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Preparation

Human recombinant monoclonal antibody purified by StrepTactin affinity chromatography from *E. coli* extracts

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

Immunogen

Cytokeratin 18

External Database Links	UniProt: P05783 Related reagents Entrez Gene: 3875 KRT18 Related reagents
Synonyms	CYK18
Specificity	<p>Human anti Human cytokeratin 18 antibody recognizes cytokeratin 18, also known as K18, CK-18, cell proliferation-inducing gene 46 protein, cell proliferation-inducing protein 46, keratin type I cytoskeletal 18 or simply keratin-18.</p> <p>Encoded by the KRT18 gene, cytokeratin 18 together with its filament partner cytokeratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in KRT18 have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for KRT18 (provided by RefSeq, Jul 2008).</p> <p>Human anti Human cytokeratin 18 antibody detects a band of 46 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti cytokeratin 18 detects a band of approximately 46 kDa in MCF7 cell lysates
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. PrecisionAb is a trademark of Bio-Rad Laboratories. His-tag is a registered trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Human IgG F(ab')₂ (STAR126...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M371797:200612'

