

Datasheet: VMA00459

**BATCH NUMBER 160712**

<b>Description:</b>	MOUSE ANTI RPS6KA3
<b>Specificity:</b>	RPS6KA3
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Preparation

Mouse monoclonal antibody isolated from ascites by saturated ammonium sulfate precipitation followed by dialysis against PBS

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)

<b>Immunogen</b>	Purified His-tagged fragment of human RPS6KA3
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P51812</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">6197</a>    RPS6KA3    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	ISPK1, MAPKAPK1B, RSK2
<b>Specificity</b>	<p><b>Mouse anti Human RPS6KA3 antibody</b> recognizes the ribosomal protein S6 kinase alpha-3, also known as MAP kinase-activated protein kinase 1b, insulin-stimulated protein kinase 1, p90-RSK 3 or ribosomal S6 kinase 2.</p> <p>RPS6KA3 encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 non-identical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signalling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Mutations in RPS6KA3 gene have been associated with Coffin-Lowry syndrome (<a href="#">CLS</a>) (provided by RefSeq, Jul 2008), additionally mutations in the RPS6KA3 gene have been associated with a non-syndromic mental retardation X-linked-19 (<a href="#">MRX19</a>).</p> <p>Mouse anti Human RPS6KA3 antibody detects a band of 83 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti RPS6KA3 detects a band of approximately 83 kDa in A431 cell lysate
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/VMA00459">https://www.bio-rad-antibodies.com/SDS/VMA00459</a> Antibody (10040)
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

**North & South America**    Tel: +1 800 265 7376  
 Fax: +1 919 878 3751

**Worldwide**    Tel: +44 (0)1865 852 700  
 Fax: +44 (0)1865 852 739

**Europe**    Tel: +49 (0) 89 8090 95 21  
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M370207:200529'

**Printed on 18 Jan 2024**

---