

Datasheet: VMA00502

**BATCH NUMBER 172192**

<b>Description:</b>	MOUSE ANTI CD324
<b>Specificity:</b>	CD324
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	AB04/4G6
<b>Isotype:</b>	IgG2b
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000
Immunofluorescence	▪			

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Learn about [how we validate our PrecisionAb range](#).** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from ascites.
<b>Buffer Solution</b>	Phosphate buffered saline.
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> ).
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml.
<b>Immunogen</b>	<i>E. coli</i> derived recombinant protein corresponding to amino acids 338-621 of human CD324

<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P12830</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b> <a href="#">999</a> CDH1 <a href="#">Related reagents</a>
<b>Synonyms</b>	CDHE, UVO
<b>Specificity</b>	<p><b>Mouse anti Human CD324 antibody</b> recognizes CD324, also known as cadherin-1, CAM 120/80, E-cadherin, epithelial cadherin or uvomorulin.</p> <p>The CDH1 gene is a classical cadherin from the cadherin superfamily. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein comprised of five extracellular cadherin repeats, a transmembrane region and a highly conserved cytoplasmic tail. Mutations in CDH1 gene are correlated with gastric, breast, colorectal, thyroid and ovarian cancer. Loss of function is thought to contribute to progression in cancer by increasing proliferation, invasion, and/or metastasis. The ectodomain of this protein mediates bacterial adhesion to mammalian cells and the cytoplasmic domain is required for internalization. Identified transcript variants arise from mutation at consensus splice sites (provided by RefSeq, Jul 2008).</p> <p>Mouse anti Human CD324 antibody detects a band of 120 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti CD324 detects a band of approximately 120 kDa in LNCaP cell lysate.
<b>Storage</b>	This product is shipped at ambient temperature. Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
<b>Guarantee</b>	12 months from date of despatch.
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/VMA00502">https://www.bio-rad-antibodies.com/SDS/VMA00502</a>
<b>Regulatory</b>	For research purposes only.

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

### Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL \(MCA691\)](#)

Product inquiries: [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M443581:250709'

**Printed on 13 Dec 2025**

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)