

## Datasheet: VMA00511

<b>Description:</b>	MOUSE ANTI eEF1A1
<b>Specificity:</b>	eEF1A1
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	F02/1E3
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse, Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid.

### Preparation

Mouse monoclonal antibody purified by affinity chromatography on Protein G from tissue culture supernatant.

### Buffer Solution

Phosphate buffered saline.

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>).

<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml.
<b>Immunogen</b>	<i>E. coli</i> derived recombinant protein corresponding to amino acids 1-462 of human eEF1A1
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P68104</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">1915</a>    eEF1A1    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	EEF1A, EF1A, LENG7
<b>Specificity</b>	<p><b>Mouse anti Human eEF1A1 antibody, clone F02/1E3</b> recognizes the elongation factor 1-alpha 1, also known as CTCL tumor antigen, cervical cancer suppressor 3, elongation factor 1 alpha subunit, elongation factor Tu, eukaryotic elongation factor 1 A-1, eukaryotic translation elongation factor 1 alpha 1-like 14, glucocorticoid receptor AF-1 specific elongation factor or leukocyte receptor cluster member 7.</p> <p>The EEF1A1 gene encodes an isoform of the alpha subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This isoform (alpha 1) is expressed in brain, placenta, lung, liver, kidney, and pancreas, and the other isoform (alpha 2) is expressed in brain, heart and skeletal muscle. This isoform is identified as an autoantigen in 66% of patients with Felty syndrome (<a href="#">Ditzel et al. 2000</a>). EEF1A1 has been found to have multiple copies on many chromosomes, some of which, if not all, represent different pseudogenes. (provided by RefSeq, Jul 2008).</p> <p>Mouse anti Human eEF1A1 antibody, clone F02/1E3 detects a band of 50 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti eEF1A1 detects a band of approximately 50 kDa in HEK293 cell lysate.
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
<b>Guarantee</b>	12 months from date of despatch.
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: Antibody (10040): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</a>
<b>Regulatory</b>	For research purposes only.

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

## Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M398764:220623'

**Printed on 12 Sep 2022**

---

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)