

Datasheet: VMA00512

Description:	MOUSE ANTI RUNX1
Specificity:	RUNX1
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody purified by affinity chromatography on Protein G from tissue culture supernatant.
Buffer Solution	Phosphate buffered saline.
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃).
Immunogen	<i>E. coli</i> derived recombinant protein corresponding to amino acids 190-453 of human RUNX1
External Database Links	UniProt: Q01196 Related reagents

Entrez Gene:[861](#) RUNX1 [Related reagents](#)

Synonyms AML1, CBFA2

Specificity **Mouse anti Human RUNX1 antibody** recognizes RUNX1, also known as AML1-EVI-1 fusion protein, Runt-related transcription factor 1, SL3/AKV core-binding factor alpha B subunit, acute myeloid leukemia 1 protein, oncogene AML-1 or polyomavirus enhancer-binding protein 2 alpha B subunit.

Core binding factor (CBF) is a heterodimeric transcription factor that binds to the core element of many enhancers and promoters. The protein encoded by RUNX1 gene represents the alpha subunit of CBF and is thought to be involved in the development of normal haematopoiesis. Chromosomal translocations involving RUNX1 are well-documented and have been associated with several types of leukemia ([Aydin et al. 2016](#)). Three transcript variants encoding different isoforms have been found for RUNX1 (provided by RefSeq, Jul 2008).

Mouse anti Human RUNX1 antibody detects a band of 55 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting Anti RUNX1 detects a band of approximately 55 kDa in MOLT-4 cell lysate.

Storage Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee 12 months from date of despatch.

Acknowledgements PrecisionAb is a trademark of Bio-Rad Laboratories.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/VMA00512>
Antibody (10040)

Regulatory For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M398767:220623'

