

## Datasheet: VMA00568

<b>Description:</b>	MOUSE ANTI RAD18
<b>Specificity:</b>	RAD18
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	AB02/3C6
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Mouse monoclonal antibody purified by affinity chromatography on Protein G from tissue culture supernatant.
<b>Buffer Solution</b>	Phosphate buffered saline.
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> ).
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	<i>E. coli</i> -derived recombinant human RAD18 (amino acids 1 - 495)

**External Database  
Links**

**UniProt:**

[Q9NS91](#) [Related reagents](#)

**Entrez Gene:**

[56852](#) RAD18 [Related reagents](#)

---

**Synonyms**

RNF73

---

**Specificity**

**Mouse anti Human RAD18 antibody** recognizes the E3 ubiquitin-protein ligase RAD18, also known as RAD18 homolog, RING finger protein 73 or postreplication repair protein RAD18.

The protein encoded by RAD18 is highly similar to *S. cerevisiae* DNA damage repair protein Rad18. Yeast Rad18 functions through its interaction with Rad6, which is an ubiquitin-conjugating enzyme required for post-replication repair of damaged DNA. Similar to its yeast counterpart, this protein is able to interact with the human homolog of yeast Rad6 protein through a conserved ring-finger motif. Mutation of this motif results in defective replication of UV-damaged DNA and hypersensitivity to multiple mutagens (provided by RefSeq, Jul 2008).

Mouse anti Human RAD18 antibody detects a bands of 80 and 90 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**

Anti RAD18 detects bands of approximately 80 and 90 kDa in K562 cell lysates.

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee**

12 months from date of despatch.

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories.

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #10040 available at:  
<https://www.bio-rad-antibodies.com/SDS/VMA00568>  
Antibody (10040)

---

**Regulatory**

For research purposes only.

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

**North & South  
America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M429508:240410'

**Printed on 10 Apr 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)