

Datasheet: VMA00679

BATCH NUMBER 100004127

Description:	MOUSE ANTI PKA C-ALPHA
Specificity:	PKA C-ALPHA
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	AB01/1D9
Isotype:	IgG1
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunoprecipitation	▪			
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - Liquid

Preparation

Mouse monoclonal antibody affinity purified on Protein G from tissue culture supernatant

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	<i>E. coli</i> -derived recombinant protein of amino acids 1-351 of human PKA C-alpha
External Database Links	<p>UniProt: P17612 Related reagents</p> <p>Entrez Gene: 5566 PRKACA Related reagents</p>
Synonyms	PKACA
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line
Specificity	<p>Mouse anti Human PKA C-alpha antibody recognizes PKA C-alpha, also known as cAMP-dependent protein kinase catalytic subunit alpha, PRKACA, PPNAD4 and protein kinase A catalytic subunit.</p> <p>Encoded by PRKACA, the ubiquitously-expressed protein PKA C-alpha is the catalytic subunit C alpha of the cyclic-AMP (cAMP)-dependent protein kinase, otherwise known as PKA (Berthon et al. 2015). PKA is a serine/threonine kinase comprised of a regulatory subunit dimer and two catalytic subunits. Activation of PKA occurs when two cAMP molecules bind to each regulatory subunit, leading to release of the catalytic subunits. The catalytic subunits act as kinases propagate downstream signaling. Through this activity, PKA is responsible for phosphorylating a wide range of substrates, and is considered an essential regulator of many cell signaling events (Turnham and Scott 2016). Somatic PRKACA mutations have been identified in tumors of the adrenal cortex (Berthon et al. 2015). PRKACA may contribute to cancer progression, with extracellular PKA catalytic subunit found in the serum of patients with cancer, and increased PRKACA transcripts observed in breast cancer (Turnham and Scott 2016).</p>
Western Blotting	Mouse anti PKA C-alpha detects a band of approximately 38 kDa in HeLa cell lysates
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles
Guarantee	12 months from date of despatch
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/VMA00679 Antibody (10040)

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M386900:210601'

Printed on 10 Apr 2024
