

Datasheet: VMA00721

**BATCH NUMBER 170727**

<b>Description:</b>	MOUSE ANTI BECLIN 1
<b>Specificity:</b>	BECLIN 1
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	OTI3F3
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Rat <b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
<b>Product Form</b>	Purified IgG - Liquid
<b>Preparation</b>	Mouse monoclonal antibody affinity purified on Protein G from ascites
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide

<b>Stabilisers</b>	1% Bovine Serum Albumin
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	HEK293-derived full length recombinant protein of human beclin 1
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">Q14457</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">8678</a>    BECN1    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	GT197
<b>Specificity</b>	<p><b>Mouse anti Human beclin 1 antibody</b> recognizes beclin 1, also known as ATG6 autophagy related 6 homolog and coiled-coil myosin-like Bcl2-interacting protein.</p> <p>Beclin 1 participates in the regulation of autophagy and has an important role in development, tumorigenesis and neurodegeneration (Zhong <i>et al.</i>, 2009 (PubMed 19270693)) (supplied by OMIM, Jul 2010).</p> <p>Mouse anti Human beclin 1 antibody detects a band of 60 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Mouse anti beclin 1 detects a band of approximately 60 kDa in Jurkat cell lysates
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10041 available at: <a href="https://www.bio-rad-antibodies.com/SDS/VMA00721">https://www.bio-rad-antibodies.com/SDS/VMA00721</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M406369:220922'

