

Datasheet: VMA00733

BATCH NUMBER 100004125

Description:	MOUSE ANTI HISTONE DEACETYLASE 3
Specificity:	HISTONE DEACETYLASE 3
Format:	Purified
Product Type:	PrecisionAb Monoclonal
Clone:	E04/3H10
Isotype:	IgG2b
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000 - 1/2000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
Product Form	Purified IgG - Liquid
Preparation	Mouse monoclonal antibody affinity purified on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide

Stabilisers

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	<i>E. coli</i> --derived recombinant protein of amino acids 1-428 of human histone deacetylase 3
External Database Links	UniProt: O15379 Related reagents Entrez Gene: 8841 HDAC3 Related reagents
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line
Specificity	Mouse anti Human histone deacetylase 3 antibody recognizes histone deacetylase 3 (HDAC3), also known as SMAP45. HDAC3 is a ubiquitously-expressed histone deacetylase, mostly located in the nucleus within transcriptional coregulator complexes (Song et al. 2019). By inhibiting the activity of multiple proteins, HDAC3 plays important roles in physiological homeostasis in many tissues (Kwon et al. 2018). In addition to its enzymatic role in catalysis of deacetylation reactions, HDAC3 can regulate chromatin accessibility and corepressor complex dynamics by acting as a chaperone or a scaffold anchorage protein, instead of an enzyme (Song et al. 2019). HDAC3 contributes to inflammation, cancer, and degenerative neurological diseases, and inhibitors of HDAC3 have been suggested as potential treatments (Cao et al. 2018). HDAC3 overexpression is associated with poor prognosis in a range of cancer types, however, the role of HDAC3 in cancer is tissue-specific. HDAC3 can also regulate responses to anticancer drugs (Kwon et al. 2018).
Western Blotting	Mouse anti histone deacetylase 3 detects a band of approximately 47 kDa in HeLa cell lysates
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles
Guarantee	12 months from date of despatch
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/VMA00733 Antibody (10040)
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M370410:200529'

Printed on 13 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)