

Datasheet: VMA00740

**BATCH NUMBER 170727**

<b>Description:</b>	MOUSE ANTI C-JUN
<b>Specificity:</b>	C-JUN
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	OT12A9
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - Liquid
<b>Preparation</b>	Mouse monoclonal antibody affinity purified on Protein A from ascites
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide
<b>Stabilisers</b>	1% Bovine Serum Albumin
<b>Immunogen</b>	<i>E. coli</i> -derived full length recombinant protein of human c-Jun
<b>External Database Links</b>	<b>UniProt:</b>

**Entrez Gene:**

---

**Specificity**      **Mouse anti Human c-Jun antibody** recognizes c-Jun, also known as Jun activation domain binding protein, activator protein 1, enhancer-binding protein AP1, proto-oncogene c-Jun and transcription factor AP-1.

JUN is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. JUN is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies (provided by RefSeq, Jul 2008).

Mouse anti Human c-Jun antibody detects bands of 43 and 48 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**      Mouse anti c-Jun detects bands of approximately 43 & 48 kDa in HeLa cell lysates

---

**Storage**      Store undiluted at -20°C, avoiding repeated freeze thaw cycles

---

**Guarantee**      12 months from date of despatch

---

**Acknowledgements**      PrecisionAb is a trademark of Bio-Rad Laboratories.

---

**Health And Safety Information**      Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/VMA00740>  
Antibody (10041)

---

**Regulatory**      For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

**North & South**      Tel: +1 800 265 7376

**America**      Fax: +1 919 878 3751

    Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

    Tel: +44 (0)1865 852 700

    Fax: +44 (0)1865 852 739

    Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

    Tel: +49 (0) 89 8090 95 21

    Fax: +49 (0) 89 8090 95 50

    Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M370414:200529'

Printed on 18 Jan 2024