

Datasheet: VMA00886

**BATCH NUMBER 100004571**

<b>Description:</b>	MOUSE ANTI CCT EPSILON
<b>Specificity:</b>	CCT EPSILON
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	AB03/2D7
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Mouse, Rat <b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
<b>Product Form</b>	Purified IgG - Liquid
<b>Preparation</b>	Mouse monoclonal antibody affinity purified on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide

## Stabilisers

---

**Approx. Protein Concentrations** IgG concentration 1.0 mg/ml

---

**Immunogen** *E. coli*-derived recombinant protein of amino acids 1-541 of human CCT epsilon

---

## External Database Links

**UniProt:**  
[P48643](#)    [Related reagents](#)

**Entrez Gene:**  
[22948](#) CCT5    [Related reagents](#)

---

**Synonyms** CCTE, KIAA0098

---

**Fusion Partners** Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line

---

**Specificity** **Mouse anti CCT epsilon antibody** recognizes T-complex protein 1 subunit epsilon, also known as CCT-epsilon and CCTE.

CCT epsilon is a component of the chaperonin-containing T-complex (TRiC), also called CCT, a large double-ring complex that appears to have the primary function of participating in protein folding ([Freund et al. 2014](#)). Each TriC subunit, named CCT1-8, is found at a fixed position within the ring. In eukaryotes, TRiC supports folding of ~10% of the cytosolic proteome, and is required for folding of the cytoskeletal protein actin ([Balchin et al. 2018](#)). Many products of genes that are deregulated in cancers are client proteins of TRiC, which suggests a role of TRiC in cancer development ([Carr et al. 2017](#)). A heritable mutation in CCT epsilon has been found to cause profound neuropathy in humans ([Spigolon et al. 2017](#)).

---

**Western Blotting** Mouse anti CCT epsilon antibody detects a band of approximately 58 kDa in HEK293 cell lysates

---

**Storage** Store undiluted at -20°C, avoiding repeated freeze thaw cycles

---

**Guarantee** 12 months from date of despatch

---

**Acknowledgements** PrecisionAb is a trademark of Bio-Rad Laboratories

---

**Health And Safety Information** Material Safety Datasheet documentation #10040 available at: Antibody (10040): <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M375107:201216'

**Printed on 09 Feb 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)