

Datasheet: VMA00969

**BATCH NUMBER 100005718**

<b>Description:</b>	MOUSE ANTI RPL6
<b>Specificity:</b>	RPL6
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Monoclonal
<b>Clone:</b>	AB01/1E12
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Mouse, Rat <b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
<b>Product Form</b>	Purified IgG - Liquid
<b>Preparation</b>	Mouse monoclonal antibody affinity purified on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide

## Stabilisers

---

**Approx. Protein Concentrations** IgG concentration 1.0 mg/ml

---

**Immunogen** *E. coli*-derived recombinant protein of amino acids 1-340 of human RPL6

---

## External Database Links

**UniProt:**

[Q02878](#)    [Related reagents](#)

**Entrez Gene:**

[6128](#) RPL6    [Related reagents](#)

---

**Synonyms** TXREB1

---

**Fusion Partners** Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line

---

**Specificity** **Mouse anti RPL6 antibody** recognizes 60S ribosomal protein L6, also known as TXREB1. RPL6 is a ribosomal protein which interacts with histone H2A and is recruited to DNA damage sites under the control of PARP. Here, RPL6 participates in the DNA damage response ([Yang et al. 2019](#)). RPL6 has been found to be upregulated in multidrug-resistant cancer cells, and could protect gastric cancer from drug-mediated apoptosis. It has been suggested that down-regulation of RPL6 in gastric cancer treatment could suppress cell growth and cell cycle progression ([Wu et al. 2011](#)). RPL6 has also been found to promote proliferation and tumorigenesis of gallbladder cancer cells through interaction with NSUN2 ([Gao et al. 2019](#)).

---

**Western Blotting** Mouse anti RPL6 detects a band of approximately 33 kDa in Jurkat cell lysates

---

**Storage** Store undiluted at -20°C, avoiding repeated freeze thaw cycles

---

**Guarantee** 12 months from date of despatch

---

**Acknowledgements** PrecisionAb is a trademark of Bio-Rad Laboratories

---

**Health And Safety Information** Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/VMA00969>  
Antibody (10040)

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

'M390065:210823'

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)