

Datasheet: VPA00011KT

Description:	ANNEXIN II ANTIBODY WITH CONTROL LYSATE
Specificity:	ANNEXIN II
Other names:	Lipocortin II
Format:	Purified
Product Type:	PrecisionAb™ Polyclonal
Isotype:	Polyclonal IgG
Quantity:	2 Westerns

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	20µl Rabbit polyclonal antibody purified by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	6x His-tagged full-length recombinant human ANXA2
External Database Links	UniProt: P07355 Related reagents Entrez Gene: 302 ANXA2 Related reagents

Synonyms	ANX2, ANX2L4, CAL1H, LPC2D
Specificity	<p>Rabbit anti Human annexin II antibody recognizes human annexin II also known as annexin A2 (ANXA2), lipocortin 2, chromobindin-8 or placental anticoagulant protein IV. Annexin II is a 339 amino acid secreted protein containing 4 annexin repeats. Annexin II is a pleiotropic calcium-dependent phospholipid-binding protein and member of the annexin family, abundantly expressed in a variety of tissues including the brain, spleen, kidney and placenta. The protein can exist in a monomer, heterodimer or heterotetramer form.</p> <p>Annexin II acts as an autocrine factor, enhancing the formation of osteoclasts and resorption of bone, and is implicated in both exocytosis and endocytosis. Annexin II is a substrate for both tyrosine and serine/threonine kinases, and is phosphorylated following growth factor stimulation and cellular activation by viral transformation (Brambilla et al. 1991).</p> <p>Rabbit anti Human annexin II antibody recognizes annexin II as a single band of ~39 kDa in A549 cell line lysates by western blotting under reducing conditions. It also cross reacts with rodent cell line lysates tested.</p>
Western Blotting	Anti annexin II antibody detects a band of approximately 39 kDa in A549 cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Lysate Composition	400µg A549 lysate lyophilized in RIPA buffer.
Lysate Reconstitution	<p>- If using DDT reconstitute the lyophilized lysate with 190µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.</p> <p>- If using BME reconstitute the lyophilized lysate with 180µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.</p> <p>Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.</p>
Storage	<p>Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.</p> <p>Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.</p>
Guarantee	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 #10561 available at:</p> <p>Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</p> <p>Lysate Material (10561): https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com 'M345712:190125'	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

Printed on 15 Mar 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)