

Datasheet: VPA00015

BATCH NUMBER 090714

Description:	SHEEP ANTI NFkB p65
Specificity:	NFkB p65
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Sheep polyclonal antibody purified by affinity chromatography

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)
0.25% Bovine Serum Albumin

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Human NFkB p65, amino acids 12-317, expressed in <i>E. coli</i>
External Database Links	<p>UniProt: Q04206 Related reagents</p> <p>Entrez Gene: 5970 RELA Related reagents</p>
Synonyms	NFKB3
Specificity	<p>Sheep anti Human NFkB p65 antibody recognizes human full length human p65 (NFkB), also known as transcription factor p65 or nuclear factor of kappa light polypeptide gene enhancer in B-cells 3. Human NFkB p65 is a 551 amino acid pleiotropic transcription factor containing a single RHD domain. Sheep anti Human NFkB p65 antibody will also detect proteolytic fragments of NFkB p65 containing the DNA binding RHD domain located in the amino-terminal region. Multiple isoforms of NFkB p65 are generated by alternative splicing. Sheep anti Human NFkB p65 antibody binds to isoforms 1 and 2 containing the intact RHD domain, but it is unclear whether it will bind to isoforms 3 and 4 which have partial differing deletions within the RHD domain.</p> <p>Sheep anti Human NFkB p65 antibody recognizes human NFkB p65 as a single band of ~65 kDa in HEK293 cell line lysates by western blotting under reducing conditions. It also cross reacts with rodent cell line lysates tested.</p>
Western Blotting	Anti NFkB p65 antibody detects a band of approximately 65 kDa in HEK293 cell lysates.
References	<ol style="list-style-type: none"> Serban, A.I. <i>et al.</i> (2015) AGEs-Induced IL-6 Synthesis Precedes RAGE Up-Regulation in HEK 293 Cells: An Alternative Inflammatory Mechanism? Int J Mol Sci. 16 (9): 20100-17. Serban, A.I. <i>et al.</i> (2016) RAGE and TGF-β1 Cross-Talk Regulate Extracellular Matrix Turnover and Cytokine Synthesis in AGEs Exposed Fibroblast Cells. PLoS One. 11 (3): e0152376.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/VPA00015 Antibody (10041)
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2504...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M370441:200529'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)