

## Datasheet: VPA00051

<b>Description:</b>	GOAT ANTI TRAF2
<b>Specificity:</b>	TRAF2
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

#### Target Species

Human

#### Species Cross Reactivity

Reacts with: Mouse, Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid

#### Preparation

Goat polyclonal antibody purified by affinity chromatography.

#### Buffer Solution

TRIS buffered saline.

#### Preservative Stabilisers

0.02% Sodium Azide (NaN<sub>3</sub>)  
0.5% Bovine Serum Albumin

#### Immunogen

Synthetic peptide sequence C-KMEAKNSYVRDD from the C-terminal region of TRAF2.

**External Database  
Links**

**UniProt:**

[Q12933](#)   [Related reagents](#)

**Entrez Gene:**

[7186](#)   TRAF2   [Related reagents](#)

---

**Synonyms**

TRAP3

---

**Specificity**

**Goat anti Human TRAF2 antibody** recognizes human tumor necrosis factor (TNF) receptor associated factor 2 (TRAF2), a member of the TNF receptor associated factor (TRAF) protein family. TRAF2 is an adapter protein and signal transducer, which acts downstream of the TNF receptor. TRAF2 directly binds with TNF receptors and other TRAF proteins, and is essential for TNF alpha mediated activation of MAPK8/JNK and NF-kappaB. TRAF2 plays a role in apoptosis, through the formation of a heterodimeric complex with TRAF1, and interaction with the IAP family members cIAP1 and cIAP2, functioning as a mediator of the anti-apoptotic signals from TNF receptors.

Goat anti Human TRAF2 recognizes TRAF2 as a band of 46 kDa in a range of human cell line lysates. Cross reactivity with rodent cell lines is also observed.

---

**Western Blotting**

Goat anti TRAF2 detects a band of approximately 46 kDa HEK293 cell lysates.

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee**

12 months from date of despatch.

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories.

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #10058 available at:  
<https://www.bio-rad-antibodies.com/SDS/VPA00051>  
Antibody (10058)

---

**Regulatory**

For research purposes only.

---

## Related Products

### Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

**North & South  
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M400172:220707'

**Printed on 25 Mar 2023**