

Datasheet: VPA00059

Description:	GOAT ANTI GRP94
Specificity:	GRP94
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Goat polyclonal antibody purified by affinity chromatography.

Buffer Solution

TRIS buffered saline.

Preservative Stabilisers

0.02% Sodium Azide (NaN₃)
0.5% Bovine Serum Albumin

Approx. Protein Concentrations

IgG concentration 0.5 mg/ml

Immunogen	Peptide with sequence C-KEGVKFDSEKTK from the internal region of the protein sequence.
External Database Links	<p>UniProt: P14625 Related reagents</p> <p>Entrez Gene: 7184 HSP90B1 Related reagents</p>
Synonyms	GRP94, TRA1
Specificity	<p>Goat anti Human GRP94 antibody recognizes GRP94 a 94 kDa glucose regulated protein (GRP94), also known as tumor rejection antigen 1 (TRA1), heat shock protein 90 kDa beta member 1 (HSP90B1) and endoplasmic reticulum resident molecular chaperone, which is a paralog of cytoplasmic hsp90. GRP94 assists in the folding of membrane-bound proteins including the toll-like receptors, certain subclasses of integrins and some secreted proteins such as IgG.</p> <p>Goat anti Human GRP94 detects GRP94 as a band of 100 kDa in a wide range of human cell line lysates tested. Cross reactivity with rodent cell lines is also observed.</p>
Western Blotting	Goat anti GRP94 detects a band of approximately 100 kDa in MCF7 cell lysates.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch.
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10058 available at: https://www.bio-rad-antibodies.com/SDS/VPA00059 Antibody (10058)
Regulatory	For research purposes only.

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...)[HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M428860:240301'

Printed on 01 Mar 2024

