

Datasheet: VPA00070

BATCH NUMBER 190914

Description:	GOAT ANTI SUPEROXIDE DISMUTASE (Cu-Zn)
Specificity:	SUPEROXIDE DISMUTASE (Cu-Zn)
Other names:	SOD1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
Product Form	Purified IgG - liquid
Preparation	Goat polyclonal antibody purified by affinity chromatography
Buffer Solution	TRIS buffered saline
Preservative	0.02% Sodium Azide (NaN ₃)

Stabilisers	0.5% Bovine Serum Albumin
Approx. Protein Concentrations	0.5 mg/ml
Immunogen	Synthetic peptide with sequence C-SRKHGGPKDEERH from the internal region of the protein.
External Database Links	<p>UniProt: P00441 Related reagents</p> <p>Entrez Gene: 6647 SOD1 Related reagents</p>
Specificity	<p>Goat anti Human superoxide dismutase (Cu-Zn) antibody recognizes superoxide dismutase (Cu-Zn) otherwise known as SOD1. Superoxide dismutases (SOD) are a class of enzymes that catalyze the dismutation of superoxide into oxygen and hydrogen peroxide. As such, they are an important antioxidant defense in nearly all cells exposed to oxygen. There are three major families of superoxide dismutase, depending on the metal cofactor: Cu-Zn (which binds both copper and zinc), Fe and Mn types (which bind either iron or manganese), and finally the Ni type, which binds nickel.</p> <p>Goat anti Human superoxide dismutase (Cu-Zn) antibody detects SOD1 as a single band of 16 kDa in a range of human cell line lysates. Cross reactivity with rodent cell lines is also observed.</p>
Western Blotting	Goat anti superoxide dismutase (Cu-Zn) detects a band of approximately 16 kDa in Jurkat cell lysates.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10058 available at: https://www.bio-rad-antibodies.com/SDS/VPA00070 Antibody (10058)
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M370474:200529'

Printed on 01 Mar 2024