

Datasheet: VPA00110

Description:	GOAT ANTI NADPH OXIDASE 1
Specificity:	NADPH OXIDASE 1
Other names:	NOX1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Goat polyclonal antibody purified by affinity chromatography.
Buffer Solution	TRIS buffered saline.
Preservative Stabilisers	0.02% Sodium Azide (NaN ₃) 0.5 % BSA
Immunogen	Peptide with the sequence C-DKATDIVTGLKQK, from the internal region of the protein sequence.
External Database Links	UniProt: Q9Y5S8 Related reagents

Entrez Gene:

[27035](#) NOX1 [Related reagents](#)

Synonyms	MOX1, NOH1
-----------------	------------

Specificity	<p>Goat anti Human NADPH oxidase 1 antibody recognizes the NADPH oxidase 1, also known as NOX1, mitogenic oxidase 1, MOX1, NADH/NADPH mitogenic oxidase subunit P65-MOX or NOH1.</p> <p>Encoded by the NOX1 gene, NADPH oxidase 1 is a member of the NADPH oxidase family of enzymes responsible for the catalytic one-electron transfer of oxygen to generate superoxide or hydrogen peroxide. Alternatively spliced transcript variants encoding multiple isoforms have been observed for NOX1 (provided by RefSeq, Nov 2012).</p> <p>Goat anti Human NADPH Oxidase 1 antibody detects a band of 72 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
--------------------	--

Western Blotting	Anti NOX1 detects a band of approximately 72 kDa in K562 cell lysates.
-------------------------	--

References	1. Waghela, B.N. <i>et al.</i> (2021) Upregulation of NOX-2 and Nrf-2 Promotes 5-Fluorouracil Resistance of Human Colon Carcinoma (HCT-116) Cells. Biochemistry (Mosc). 86 (3): 262-74.
-------------------	---

Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
----------------	---

Guarantee	12 months from date of despatch.
------------------	----------------------------------

Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
-------------------------	---

Health And Safety Information	Material Safety Datasheet documentation #10058 available at: Antibody (10058): https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf
--------------------------------------	--

Regulatory	For research purposes only.
-------------------	-----------------------------

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M400407:220711'

Printed on 11 Jul 2022