

Datasheet: VPA00200

BATCH NUMBER 150615

Description:	SHEEP ANTI ALPHA 2 MACROGLOBULIN
Specificity:	ALPHA 2 MACROGLOBULIN
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunoprecipitation	▪			
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Sheep polyclonal antibody purified by affinity chromatography on Protein G.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃).

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Purified human alpha 2 macroglobulin
External Database Links	<p>UniProt: P01023 Related reagents</p> <p>Entrez Gene: 2 A2M Related reagents</p>
Synonyms	CPAMD5
Specificity	<p>Sheep anti Human alpha 2 macroglobulin antibody recognizes the alpha 2 macroglobulin, also known as C3 and PZP-like alpha-2-macroglobulin domain-containing protein 5 and alpha-2-M.</p> <p>Alpha 2 macroglobulin is a protease inhibitor and cytokine transporter. It inhibits many proteases, including trypsin, thrombin and collagenase. Alpha 2 macroglobulin is implicated in Alzheimer disease (AD) due to its ability to mediate the clearance and degradation of A-beta, the major component of beta-amyloid deposits (provided by RefSeq, Jul 2008).</p> <p>Sheep anti Human alpha 2 macroglobulin antibody detects a band of 180 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti alpha 2 macroglobulin detects a band of approximately 180 kDa in HepG2 cell lysates.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch.
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/VPA00200 Antibody (10040)
Regulatory	For research purposes only.

Related Products

Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2504...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M404432:220824'

Printed on 01 Mar 2024
