

Datasheet: VPA00215

BATCH NUMBER 151204

Description:	RABBIT ANTI EIF3A
Specificity:	EIF3A
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Rabbit polyclonal antibody purified by affinity chromatography

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

Immunogen KLH conjugated synthetic peptide between 1219-1247 amino acids from the C-terminal region of human EIF3A

External Database Links

UniProt:

[Q14152](#) [Related reagents](#)

Entrez Gene:

[8661](#) EIF3A [Related reagents](#)

Synonyms EIF3S10, KIAA0139

Specificity **Rabbit anti Human EIF3A antibody** recognizes the eukaryotic translation initiation factor 3 subunit A, also known as EIF3, p180 subunit, centrosomin homolog, cytoplasmic protein p167, eIF-3-theta, eIF3 p167, eIF3 p180, eIF3 p185 and eukaryotic translation initiation factor 3 subunit 10.

Rabbit anti Human Eukaryotic translation initiation factor 3 subunit A antibody detects a band of 166 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting Anti EIF3A detects a band of approximately 166 kDa in PC-12 cell lysates

Storage Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee 12 months from date of despatch

Acknowledgements PrecisionAb is a trademark of Bio-Rad Laboratories.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/VPA00215>
Antibody (10040)

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M370579:200529'

Printed on 13 Aug 2023
