

Datasheet: VPA00276

**BATCH NUMBER 160415**

<b>Description:</b>	RABBIT ANTI Hsc70
<b>Specificity:</b>	Hsc70
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	<p>Reacts with: Mouse</p> <p><b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
<b>Product Form</b>	Ig fraction - liquid
<b>Preparation</b>	Rabbit Ig fraction prepared by ammonium sulphate precipitation
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )

**Immunogen** KLH-conjugated synthetic peptide corresponding to aa 82-110 of human Hsc70

**External Database**

**Links**

**UniProt:**

[P11142](#)

[Related reagents](#)

**Entrez Gene:**

[3312](#)

HSPA8

[Related reagents](#)

**Synonyms**

HSC70, HSP73, HSPA10

**Specificity**

**Rabbit anti Human Hsc70 antibody** recognizes Hsc70 also known as HSPA8, HSC71, N-myristoyltransferase inhibitor protein 71, constitutive heat shock protein 70, epididymis luminal protein 33, epididymis secretory sperm binding protein Li 72p, heat shock 70kd protein 10, heat shock cognate protein 54 and lipopolysaccharide-associated protein 1.

The Hsc70 gene encodes a member of the heat shock protein 70 family, which contains both heat-inducible and constitutively expressed members. This protein belongs to the latter group, which are also referred to as heat-shock cognate proteins. It functions as a chaperone, and binds to nascent polypeptides to facilitate correct folding. It also functions as an ATPase in the disassembly of clathrin-coated vesicles during transport of membrane components through the cell. Alternatively spliced transcript variants encoding different isoforms have been found for Hsc70 (provided by RefSeq, Aug 2011).

Rabbit anti Human Hsc70 protein antibody detects a band of 75 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

**Western Blotting**

Anti Hsc70 detects a band of approximately 75 kDa in HEK293 cell lysate

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles

**Guarantee**

12 months from date of despatch

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories.

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at:

<https://www.bio-rad-antibodies.com/SDS/VPA00276>

Antibody (10040)

**Regulatory**

For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M370637:200529'

**Printed on 13 Aug 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)