

Datasheet: VPA00294

Description:	RABBIT ANTI SIRT1
Specificity:	SIRT1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Rabbit polyclonal antibody purified by affinity chromatography.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃).

Immunogen

KLH-conjugated synthetic peptide corresponding to aa 587-621 of human SIRT1

**External Database
Links**

UniProt:

[Q96EB6](#) [Related reagents](#)

Entrez Gene:

[23411](#) SIRT1 [Related reagents](#)

Synonyms

SIR2L1

Specificity

Rabbit anti Human SIRT1 antibody recognizes SIRT1, also known as NAD-dependent deacetylase sirtuin-1, SIR2-like protein 1, SIR2alpha, hSIR2, hSIRT1, regulatory protein SIR2 homolog 1 and sirtuin type 1.

The SIRT1 gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. Human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity ([Carafa et al. 2016](#)). The protein encoded by SIRT1 is included in class I of the sirtuin family. Alternative splicing results in multiple transcript variants (provided by RefSeq, Dec 2008).

Rabbit anti Human SIRT1 antibody detects a band of 120 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti SIRT1 detects a band of approximately 120 kDa in NIH/3T3 cell lysate.

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee

12 months from date of despatch.

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories.

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/VPA00294>
Antibody (10040)

Regulatory

For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M401139:220714'

Printed on 25 Mar 2023

