

Datasheet: VPA00302

BATCH NUMBER 160411

| | |
|----------------------|------------------------|
| Description: | RABBIT ANTI BMP-4 |
| Specificity: | BMP-4 |
| Format: | Purified |
| Product Type: | PrecisionAb Polyclonal |
| Isotype: | Polyclonal IgG |
| Quantity: | 100 µl |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ■ | | | 1/1000 |

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

| | |
|---------------------------------|--|
| Target Species | Human |
| Product Form | Purified IgG - liquid |
| Preparation | Rabbit polyclonal antibody purified by affinity chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 0.09% Sodium Azide (NaN ₃) 2% Sucrose |
| Immunogen | Synthetic peptide corresponding to the C-terminal region of human BMP-4 |
| External Database Links | UniProt: P12644 Related reagents |

Entrez Gene:[652](#) [BMP4](#) [Related reagents](#)

| | |
|--------------------------------------|--|
| Synonyms | BMP2B, DVR4 |
| Specificity | <p>Rabbit anti Human BMP-4 antibody recognizes bone morphogenetic protein 4 (BMP-4) also known as bone morphogenetic protein 2B (BMP-2B).</p> <p>BMP-4 is a member of the bone morphogenetic protein family which is part of the transforming growth factor-beta superfamily. Bone morphogenetic proteins were originally identified by an ability of demineralized bone extract to induce endochondral osteogenesis in vivo at an extraskeletal site. This particular family member plays an important role in the onset of endochondral bone formation in humans, and a reduction in expression has been associated with a variety of bone diseases, including the heritable disorder Fibrodysplasia Ossificans Progressiva. Alternative splicing in the 5' untranslated region of the BMP-4 gene sequence has been described and three variants reported, all encoding an identical protein (provided by RefSeq, Jul 2008).</p> <p>Rabbit anti Human BMP-4 antibody detects a band of 47 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p> |
| Western Blotting | Anti BMP-4 detects a band of approximately 47 kDa in A549 cell lysate |
| Storage | Store undiluted at -20°C, avoiding repeated freeze thaw cycles |
| Guarantee | 12 months from date of despatch |
| Acknowledgements | PrecisionAb is a trademark of Bio-Rad Laboratories. |
| Health And Safety Information | Material Safety Datasheet documentation #10045 available at: https://www.bio-rad-antibodies.com/SDS/VPA00302 Antibody (10045) |
| Regulatory | For research purposes only |

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M370664:200529'

Printed on 13 Aug 2023

