

Datasheet: VPA00328

**BATCH NUMBER 160309**

<b>Description:</b>	RABBIT ANTI PPP2R2A
<b>Specificity:</b>	PPP2R2A
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	<p>Reacts with: Mouse</p> <p><b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )

**Immunogen** KLH-conjugated synthetic peptide corresponding to aa 126-160 of human PPP2R2A

---

**External Database**

**Links**

**UniProt:**

[P63151](#) [Related reagents](#)

**Entrez Gene:**

[5520](#) PPP2R2A [Related reagents](#)

---

**Specificity**

**Rabbit anti Human PPP2R2A antibody** recognizes PPP2R2A, also known as protein phosphatase 2 (formerly 2A), regulatory subunit B (PR52) alpha isoform, 55 kDa regulatory subunit B, alpha isoform and serine/threonine-protein phosphatase 2A.

The product of the PPP2R2A gene belongs to the phosphatase 2 regulatory subunit B family. Protein phosphatase 2 is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, which associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. The PPP2R2A gene encodes an alpha isoform of the regulatory subunit B55 subfamily. Alternatively spliced transcript variants have been described (provided by RefSeq, Apr 2010).

Rabbit anti Human PPP2R2A antibody detects a band of 52 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**

Anti PPP2R2A detects a band of approximately 52 kDa in HeLa cell lysate

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories.

---

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at:  
Antibody (10040): <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M370683:200529'

