

Datasheet: VPA00346

BATCH NUMBER 160309

Description:	RABBIT ANTI UBE4B
Specificity:	UBE4B
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Ig fraction - liquid
Preparation	Rabbit Ig fraction prepared by ammonium sulphate precipitation
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	KLH-conjugated synthetic peptide corresponding to aa 112-142 of human UBE4B
External Database Links	UniProt: O95155 Related reagents

Entrez Gene:[10277](#) UBE4B [Related reagents](#)

Synonyms HDNB1, KIAA0684, UFD2

Specificity **Rabbit anti Human UBE4B antibody** recognizes UBE4B, also known as UBOX3, UFD2A-III/UBE4B-III splice isoform, homologous to yeast UFD2, homozygously deleted in neuroblastoma-1, ubiquitin conjugation factor E4B, ubiquitin-fusion degradation protein 2 and ubiquitination factor E4B (homologous to yeast UFD2).

The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. The UBE4B gene encodes an additional conjugation factor, E4, which is involved in multi-ubiquitin chain assembly. UBE4B is also the strongest neuroblastoma tumor suppressor gene. Alternatively spliced transcript variants encoding distinct isoforms have been found for UBE4B (provided by RefSeq, Jul 2008).

Rabbit anti Human UBE4B antibody detects a band of 140 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting Anti UBE4B detects a band of approximately 140 kDa in K562 cell lysate

Storage Store undiluted at -20°C, avoiding repeated freeze thaw cycles

Guarantee 12 months from date of despatch

Acknowledgements PrecisionAb is a trademark of Bio-Rad Laboratories.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/VPA00346>
Antibody (10040)

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M370699:200529'

Printed on 13 Aug 2023

