

Datasheet: VPA00388KT

| | |
|----------------------|------------------------------------|
| Description: | HBS1L ANTIBODY WITH CONTROL LYSATE |
| Specificity: | HBS1L |
| Format: | Purified |
| Product Type: | PrecisionAb Polyclonal |
| Isotype: | Polyclonal IgG |
| Quantity: | 2 Westerns |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ▪ | | | 1/1000 |

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

| | |
|---------------------------------|--|
| Target Species | Human |
| Species Cross Reactivity | Reacts with: Mouse N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. |
| Product Form | Purified IgG - liquid |
| Preparation | 20µl Rabbit polyclonal antibody purified by affinity chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin |
| Immunogen | KLH-conjugated synthetic peptide corresponding to aa 460-486 of human HBS1L |

| | |
|--------------------------------|--|
| External Database Links | UniProt: Q9Y450 Related reagents |
| | Entrez Gene: 10767 HBS1L Related reagents |
| Synonyms | HBS1, KIAA1038 |
| Specificity | <p>Rabbit anti Human HBS1L antibody recognizes HBS1L, also known as ERF3-similar protein, Hsp70 subfamily B suppressor 1-like protein and eRF3 family member.</p> <p>The HBS1L gene encodes a member of the GTP-binding elongation factor family. It is expressed in multiple tissues with the highest expression in heart and skeletal muscle. The inter-genic region of HBS1L and the MYB gene has been identified to be a quantitative trait locus (QTL) controlling fetal hemoglobin level, and this region influences erythrocyte, platelet, and monocyte counts as well as erythrocyte volume and hemoglobin content. DNA polymorphisms at this region associate with fetal hemoglobin levels and pain crises in sickle cell disease. A single nucleotide polymorphism in exon 1 of HBS1L is significantly associated with severity in beta-thalassemia/Hemoglobin E. Multiple alternatively spliced transcript variants encoding different protein isoforms have been found for HBS1L (provided by RefSeq, May 2009).</p> <p>Rabbit anti Human HBS1L antibody detects a band of 70 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p> |
| Western Blotting | Anti HBS1L detects a band of approximately 70 kDa in HEK293 cell lysate |
| Instructions For Use | Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page. |
| Lysate Composition | 400µg HEK293 lysate lyophilized in RIPA buffer |
| Lysate Reconstitution | <p>- If using DDT reconstitute the lyophilized lysate with 190µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.</p> <p>- If using BME reconstitute the lyophilized lysate with 180µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.</p> <p>Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.</p> |
| Storage | <p>Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.</p> <p>Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.</p> |
| Guarantee | As supplied, 12 months from date of despatch |
| Acknowledgements | PrecisionAb™ is a trademark of Bio-Rad Laboratories |

**Health And Safety
Information**

Material Safety Datasheet documentation #10041 #10561 available at:
Antibody (10041): <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>
Lysate Material (10561): <https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.

Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M351678:190318'

Printed on 10 Feb 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)