

## Datasheet: VPA00391

<b>Description:</b>	RABBIT ANTI CKAP4
<b>Specificity:</b>	CKAP4
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> ) 2% Sucrose
<b>Immunogen</b>	Synthetic peptide corresponding to the central region of human CKAP4
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">Q07065</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b> <a href="#">10970</a> CKAP4 <a href="#">Related reagents</a>

<b>Specificity</b>	<b>Rabbit anti Human CKAP4 antibody</b> recognizes the CKAP4, also known as cytoskeleton-associated protein 4, 63 kDa cytoskeleton-linking membrane protein or endoplasmic reticulum/Golgi intermediate compartment, type-II transmembrane protein p63.
--------------------	---

Rabbit anti Human CKAP4 antibody detects a band of 64 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

<b>Western Blotting</b>	Anti CKAP4 detects a band of approximately 64 kDa in HEK293 cell lysate
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10045 available at: Antibody (10045): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf</a>
<b>Regulatory</b>	For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M370737:200529'

**Printed on 11 Aug 2020**

---

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)