

Datasheet: VPA00397

**BATCH NUMBER 160411**

<b>Description:</b>	RABBIT ANTI ARIH1
<b>Specificity:</b>	ARIH1
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	<p>Reacts with: Mouse</p> <p><b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography on Protein A
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	<p>0.09% Sodium Azide (NaN<sub>3</sub>)</p> <p>2% Sucrose</p>

**Immunogen** Synthetic peptide corresponding to the C-terminal region of human ARIH1

**External Database**

**Links**

**UniProt:**

[Q9Y4X5](#)

[Related reagents](#)

**Entrez Gene:**

[25820](#)

ARIH1

[Related reagents](#)

**Synonyms**

ARI, MOP6, UBCH7BP

**Specificity**

**Rabbit anti Human ARIH1 antibody** recognizes ARIH1, also known as Ariadne RBR E3 ubiquitin protein ligase 1, ubiquitin-conjugating enzyme E2-binding protein 1 and UBCH7BP.

Rabbit anti Human ARIH1 antibody detects a band of 64 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

**Western Blotting**

Anti ARIH1 detects a band of approximately 64 kDa in HEK293 cell lysate

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles

**Guarantee**

12 months from date of despatch

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories

**Health And Safety Information**

Material Safety Datasheet documentation #10045 available at:  
<https://www.bio-rad-antibodies.com/SDS/VPA00397>  
Antibody (10045)

**Regulatory**

For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M370741:200529'

Printed on 17 Sep 2024