

Datasheet: VPA00397

BATCH NUMBER 172049

Description:	RABBIT ANTI ARIH1
Specificity:	ARIH1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Learn about [how we validate our PrecisionAb range](#). Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Rabbit polyclonal antibody purified by affinity chromatography on Protein A.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)
2% Sucrose

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml.
Immunogen	Synthetic peptide corresponding to the C-terminal region of human ARIH1
External Database Links	<p>UniProt: Q9Y4X5 Related reagents</p> <p>Entrez Gene: 25820 ARIH1 Related reagents</p>
Synonyms	ARI, MOP6, UBCH7BP
Specificity	<p>Rabbit anti Human ARIH1 antibody recognizes ARIH1, also known as Ariadne RBR E3 ubiquitin protein ligase 1, ubiquitin-conjugating enzyme E2-binding protein 1 and UBCH7BP.</p> <p>Rabbit anti Human ARIH1 antibody detects a band of 64 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti ARIH1 detects a band of approximately 64 kDa in HEK293 cell lysate.
Storage	This product is shipped at ambient temperature. Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch.
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10045 available at: https://www.bio-rad-antibodies.com/SDS/VPA00397
Regulatory	For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M444226:250709'

Printed on 29 Jan 2026